







Littmann®

Stethoscopes

Choose the right stethoscope to help you perform at your best

	Littmann® CORE Digital Stethoscope	Littmann® Master Cardiology™ Stethoscope	Littmann® Cardiology IV™ Stethoscope	Littmann® Classic III™ Stethoscope	Littmann® Classic II™ Pediatric Stethoscope	Littmann® Lightweight II S.E. Stethoscope
	Switch from analog listening to digital. Connects to Eko software on a smart device** to capture, store and share data. Chestpiece diameter double-sided Adult: 1.7"/4.3cm Pediatric: 1.4"/3.4cm	Best acoustic response in Littmann mechanical line. Chestpiece diameter single-sided 1.9"/4.8cm	Diagnose and make decisions with confidence. Chestpiece diameter double-sided Adult: 1.7"/4.3cm Pediatric: 1.4"/3.4cm	Monitor and assess a wide range of patients. Chestpiece diameter double-sided Adult: 1.7"/4.3cm Pediatric: 1.4"/3.4cm	Specially designed for the smallest patients. Chestpiece diameter double-sided Pediatric: 1.4"/3.4cm, 1.1"/2.8cm Infant: 1.1"/2.8cm, 0.7"/1.7cm	Teardrop-shaped chestpiece fits under blood pressure cuff. Chestpiece diameter double-sided 2.0"x1.8"/5.2cm x 4.5cm 1.4"x1.3"/3.5cm x 3.2cm
						
Acoustic performance level	10+	10	9	8	7	6
Length	27"/69cm	27"/69cm	27"/69cm 22"/56cm	27"/69cm	28"/71cm	28"/71cm
Weight	200g+	161g-200g	161g-200g	131g-160g	<100g, 100g-130g	100g-130g
Warranty	1-3 years	7 years	7 years	5 years	3 years	2 years
Diaphragm tunable	●	●	●	●		●
Up to 40x amplification*	●					
Noise cancellation	●					
Mobile connectivity	●					
Dual-lumen tubing	●	●	●			
Additional special finishes	●	●	●	●	●	
Convertible open bell	●		●	●		
Soft sealing eartips	●	●	●	●	●	●
Dual-sided chestpiece	●		●	●	●	●

*At peak frequency vs. analog mode.

**Subscription may be required.

© Solventum 2025. Solventum, the S logo, Littmann, the L Littmann logo, the L, and Cardiology IV are marks and/or registered marks of Solventum or its affiliates. Other trademarks are the property of their respective owners. 70-2013-1757-8

Learn more at www.littmann.com

 solventum