






Choose the right stethoscope to help you perform at your best

| | Littmann® Master Cardiology™ Stethoscope | Littmann® Cardiology IV™ Stethoscope | Littmann® Classic III™ Stethoscope | Littmann® Classic II™ Pediatric Stethoscope | Littmann® Lightweight II S.E. Stethoscope |
|-----------------------------|---|---|---|---|---|
| | Best acoustic response in Littmann mechanical line. Chestpiece diameter single-sided 1.9"/4.8cm | Diagnose and make decisions with confidence. Chestpiece diameter double-sided Adult: 1.7"/4.3cm Pediatric: 1.4"/3.4cm | Monitor and assess a wide range of patients. Chestpiece diameter double-sided Adult: 1.7"/4.3cm Pediatric: 1.4"/3.4cm | Specially designed for the smallest patients. Chestpiece diameter double-sided Pediatric: 1.4"/3.4cm, 1.1"/2.8cm Infant: 1.1"/2.8cm, 0.7"/1.7cm | Teardrop-shaped chestpiece fits under blood pressure cuff. Chestpiece diameter double-sided 2.0"x1.8"/5.2cm x 4.5cm 1.4"x1.3"/3.5cm x 3.2cm |
| |  |  |  |  |  |
| Acoustic performance level | 10 | 9 | 8 | 7 | 6 |
| Length | 27"/69cm | 27"/69cm 22"/56cm | 27"/69cm | 28"/71cm | 28"/71cm |
| Weight | 161g-200g | 161g-200g | 131g-160g | <100g, 100g-130g | 100g-130g |
| Warranty | 7 years | 7 years | 5 years | 3 years | 2 years |
| Diaphragm tunable | ● | ● | ● | | ● |
| Dual-lumen tubing | ● | ● | | | |
| Additional special finishes | ● | ● | ● | ● | |
| Convertible open bell | | ● | ● | | |
| Soft sealing eartips | ● | ● | ● | ● | ● |
| Dual-sided chestpiece | | ● | ● | ● | ● |

*At peak frequency vs. analog mode.
**Subscription may be required.

© Solventum 2025. Solventum, the S logo, Littmann, the L Littmann logo, the L, and Cardiology IV are marks and/or registered marks of Solventum or its affiliates. Other trademarks are the property of their respective owners. 70-2013-1757-8