

Littmann®

Stethoscopes

Choose the right stethoscope to help you perform at your best

	Littmann® Master Cardiology™ Stethoscope	Littmann® Cardiology IV™ Stethoscope	Littmann® Classic III™ Stethoscope	Littmann® Classic II™ Paediatric and Infant Stethoscopes	Littmann® Lightweight II S.E. Stethoscope
	Best acoustic response in Littmann mechanical line. Chestpiece diameter single-sided 4.8cm	Diagnose and make decisions with confidence. Chestpiece diameter double-sided Adult: 4.3cm Paediatric: 3.4cm	Monitor and assess a wide range of patients. Chestpiece diameter double-sided Adult: 4.3cm Paediatric: 3.4cm	Specially designed for the smallest patients. Chestpiece diameter double-sided Paediatric: 3.4cm, 2.8cm Infant: 2.8cm, 1.7cm	Teardrop-shaped chestpiece fits under blood pressure cuff. Chestpiece diameter double-sided 5.2cm x 4.5cm 3.5cm x 3.2cm
					
Acoustic performance level	10	9	8	7	6
Approx. length	69cm	56cm	69cm	71cm	71cm
Approx. weight	185g	177g	150g	105g/95g	118g
Warranty	7 years	7 years	5 years	3 years	2 years
Tunable diaphragm	●	●	●		●
Dual-lumen tubing	●	●			
Additional special finishes	●	●	●	●	
Convertible open bell		●	●		
Soft sealing eartips	●	●	●	●	●
Dual-sided chestpiece		●	●	●	●

© Solventum 2026. Solventum, the S logo and Littmann, the L Littmann logo, the L, Master Cardiology, Cardiology IV and Classic III are trademarks of Solventum or its affiliates. Other trademarks are the property of their respective owners. OMG2101326.

Learn more at www.littmann.com/en-za

 solventum